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Instructions for use
and safety recommendations

Intensiv PrepTwins Pat. 2564803

GB

Product description

- Diamond Polymer Finisher highly loaded in 3 different diamond grains each to refine preparation
- Form congruent instruments to the preparation diamond instruments previously applied
- Instruments in three diamond grains each, 60µm (brown), 40µm (red) and 15µm (yellow)
- The metal core below the polymer is coated with the same diamond grain as the diamond polymer mixture
- Forms: ball, cylinder, flame, football, torpedo, tapered con
- Right angle instrument, low speed.
- Maximum diameter of the working part: ISO 023
- Instruments should be stored in the original packaging.

Intended use

Finishing of the specific created preparation.

Instructions for use

- Application by dentist only.
- Apply Intensiv PrepTwins for the finishing of cavity and prosthetic preparations after using rotating diamond instruments.
- It is appropriate to choose Intensiv PrepTwins of the same shape as that of the previously used diamond instrument, in order to avoid altering the preparation.
- The choice of grit-size (60, 40 or 15µm) depends on the necessary modeling or finishing degree for cavity and prosthetic preparations.
- Insert the instrument into the chuck until safety “click” is noted
- Set the contra-angle at the required working speed up to a maximum of 20'000 rpm before applying the instrument on the area to be treated.
- Activate water spray during the entire treatment (minimum 50 ml/min) for cooling and instruments cleaning.
- Ensure that there is consistent water flow, regulate in case of water excess or shortage.
- Apply no pressure during usage.
- It is recommended to use a dental dam during treatment.
- The number of possible treatments depends on each treatment time and on the degree of finishing
- Take a new instrument when the grinding performance is no longer effective.
- Once the preparation has been completed, remove the instrument from the treated area and bring it to a standstill.
- It is recommended to wear gloves and safety glasses.

Maintenance and Sterilization

- Instruments are delivered non sterile. They must be disinfected and sterilized prior to first use on the patient and disinfected, cleaned, and sterilized after each use.
- Immerse the instruments in a disinfection container immediately after use.
- Manual cleaning and disinfection: carry out ultrasonic assisted cleaning and disinfection with cleaning agent and disinfectant suitable for polymeric polishers. Exposure time and concentration, as indicated by the manufacturer, are to be respected. Rinse the instruments afterwards under running water. Dry with single use sterile lint-free cloth or medical compressed air.
- Automatic cleaning and disinfection: carry out automatic cleaning and disinfection with thermal disinfectior, using a cleaning agent and disinfectant suitable for polymeric polishers. Cleaning and disinfection program as indicated by the manufacturer in the operating instructions. Check the instruments are dry, and if necessary dry with single use sterile lint-free cloth or medical compressed air.
- Perform visual check of all instruments for cleanliness, integrity and functionality, if necessary by using an illuminated magnifying glass (3-6 dioptres).
- If no residues are visible continue to sterilization. If residues are visible repeat the cleaning and disinfection procedure, with the aid of a soft brush.
- Dispose of instruments in the event of visible defects
- Sterilization:
 - Sterilization must be carried out according to validated procedures

- Use a single-pulsed or fractionated vacuum autoclave and sub-vacuum drying
- Sterilizable with vacuum and extra pressure steam cycles at 134°C for at least 4 min
- Do not use a chemiclave or a hot air disinfectior
- If several instruments are to be sterilized in the autoclave, ensure that the maximum weight allowed by the autoclave manufacturer is not exceeded

Risk Warnings

- Avoid jamming or levering actions when rotating, as this increases the risk of instrument breakage.
- Never exceed the specified maximum speed so as to avoid instrument breakage.
- Breakage of the polymeric coating could occur. To minimize risks comply with usage (no pressure, recommended speed, water cooling) and maintenance indications.
- Application of excessive force and/or torsion during usage may cause detachment of the polymeric coating.
- Using a disinfection solution which is not suitable for polymeric polishers may cause damage or detachment of the polymeric coating.
- Temperatures above 134°C should be avoided not to jeopardize the durability of the instrument.
- To ensure traceability of the instruments during their entire application, we recommend keeping the packaging.

Intensiv PrepTwins, Trial Kit, 6 different shapes, 15 µm, REF RA PT2415/G
Intensiv PrepTwins, Trial Kit, 6 different shapes, 40 µm, REF RA PT2440/G
Intensiv PrepTwins, Trial Kit, 6 different shapes, 60 µm, REF RA PT2460/G

Intensiv PrepTwins, 6 x ball, 60 µm, REF RA PT801/G
Intensiv PrepTwins, 6 x football, 60 µm, REF RA PT368/G
Intensiv PrepTwins, 6 x cylinder, 60 µm, REF RA PT882/G
Intensiv PrepTwins, 6 x flame, 60 µm, REF RA PT862/G
Intensiv PrepTwins, 6 x tapered cone, 60 µm, REF RA PT847KR/G
Intensiv PrepTwins, 6 x torpedo, 60 µm, REF RA PT877K/G

Intensiv PrepTwins, 6 x ball, 40 µm, REF RA PT4801/G
Intensiv PrepTwins, 6 x football, 40 µm, REF RA PT4368/G
Intensiv PrepTwins, 6 x cylinder, 40 µm, REF RA PT4882/G
Intensiv PrepTwins, 6 x flame, 40 µm, REF RA PT4862/G
Intensiv PrepTwins, 6 x tapered cone, 40 µm, REF RA PT4847KR/G
Intensiv PrepTwins, 6 x torpedo, 40 µm, REF RA PT4877K/G

Intensiv PrepTwins, 6 x ball, 15 µm, REF RA PT5801/G
Intensiv PrepTwins, 6 x football, 15 µm, REF RA PT5368/G
Intensiv PrepTwins, 6 x cylinder, 15 µm, REF RA PT5882/G
Intensiv PrepTwins, 6 x flame, 15 µm, REF RA PT5862/G
Intensiv PrepTwins, 6 x tapered cone, 15 µm, REF RA PT5847KR/G
Intensiv PrepTwins, 6 x torpedo, 15 µm, REF RA PT5877K/G

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EVE Ernst Vetter GmbH