

# Intensiv DiaTweezer

Tweezers with diamond coated tips

Good retention for gripping and holding of delicate items is often required to prevent slipping. Tweezers with open tips are useful for holding hollow objects such as crowns.

## Indications

- For gripping and securely holding small extraoral objects in dentistry e.g. a crown or an inlay/onlay

## Benefits

- Firm and secure grip through the diamond-coated inner side of the instrument tips
- Secure hold of small hollow bodies with open tweezers through fully diamond-coated instrument tips

Intensiv DiaTweezer College



REF 1001



REF 1001A

Intensiv DiaTweezer Meriam



REF 1002



REF 1002A



REF 1001: diamond-coating all around



REF 1002A: inside diamond-coating only

Classic

## Indications

- Safe gripping and holding of small extraoral objects in dentistry (e.g. a crown or an inlay/onlay)

## Benefits

- Comfortable round handle design
- Gentle grip thanks to integrated closing stop
- Firm and secure grip through the diamond-coated inner side of the tweezer tips
- Secure hold of small hollow bodies with open tweezers through fully diamond-coated tweezer tips


Intensiv DiaTweezer College



REF 1005



REF 1005A



Unique integrated closing stop



New