

# Intensiv HygienicTray


New

Resistant high quality stainless steel instrument tray for the arrangement of rotating instruments for routine treatment protocols


Instrument trays for the permanent availability of all rotating instruments for routine treatment protocols have become an indispensable support in dental clinics. The Intensiv HygienicTray offers the appropriate solution. Dental clinic processes have become standardized by the availability of the complete set of commonly used rotating instruments.



REF HT1000



REF HT2000



REF HT3000

## Product Description

- Can be disassembled into different parts made of high quality stainless steel
- Without hard to reach fissures and hollows
- Rounded edges
- Height: suitable for standard sterilization trays

### REF HT1000:

- 2 separable parts: foldable cover; instrument tray
- 13 FG and 5 RA instrument holders
- Ergonomic instrument access due to low and inclined tray position
- Dimensions: 50x40x30 mm (LxWxH)

### REF HT2000:

- 2 separable parts: foldable cover; instrument tray
- 26 FG and 10 RA instrument holders
- Ergonomic instrument access due to low and inclined tray position
- Dimensions: 104x44x30 mm (LxWxH)

### REF HT3000:

- 2 separable parts
- 36 FG and 18 RA instrument holders
- Ergonomic instrument access due to low and inclined tray position
- Dimensions: 150x44x30 mm (LxWxH)

## Application

- Organization of standard instruments for all routine treatment protocols
- Storage and transport of rotary instruments

## Benefits

- Organization and availability of commonly used rotating instruments
- Easy instrument removal due to irregular arranged position of instruments
- Hygienic due to 2 individual parts without hard to reach fissures and hollows
- Durable, resistant and sterilizable stainless steel