

Intensiv IPR-DistanceControl

Manual stainless steel instrument to support the control of tooth distances created during IPR protocols

Indications

- Support for the control of the teeth distance in the various stages of the IPR technique in orthodontics.

Benefits

- Secure application due to the handle
- Support thanks to the different thicknesses
- Large number of controls of distances, thanks to the combination of 8 instruments

Dr. Francesco Garino, Turin, Italy



Product Description

- Intensiv IPR-DistanceControl is a manual
- metal instrument
- It consists of an eight-part instrument set made of stainless steel with grip (fingergrip)
- The thickness indication is marked on each instrument
- Removable plastic rivet to secure the 8 instruments
- Height / Length of metallic element (excl. grip): 4 / 40 mm



REF IPR-DC Set