


Intensiv Metal Diamond Strips Trial Kit

Contouring, Finishing and Polishing of
aesthetical restorations in proximal surfaces

Intensiv ProxoStrip, Intensiv ProxoStrip Anterior, Intensiv ProxoContour Coars and Intensiv ProxoPolish Plus are unique diamond-coated metal strips. The length of the strips makes it easier to support the fingers during treatment. The strip height of 2,5mm and the free zone between the working areas maintain the existing contact point.



Product Description

- 4 different Metal strip, diamond-coated on one side only with two working parts
- Diamond-free zone between the diamond-coated surfaces
- Height: 2,5mm
- Thickness: 0,05 mm
- Sterilizable and reusable

Trial Kit REF PXT/6

Content:

- 2 Intensiv ProxoStrip
- 1 Intensiv ProxoStrip Anterior
- 2 Intensiv ProxoContour Coarse
- 1 Intensiv ProxoPolish