

# Intensiv Ortho-Strips® Sets

New

## Set for Interproximal Reduction

Intensiv Ortho-Strips Tray 400 Set is for the aligner treatment. It also contains the new, patented Intensiv Ortho-Strips System, Central. Thanks to the diamond coating in the central operating area, the Intensiv Ortho-Strips System, Central allows for high precision in the approximal reduction, finishing and polishing of enamel, respecting the original morphology and dental tissue.



REF OST400Set01/6



REF OST400Set02/6



REF OST400Set03/6

## Product Description

### REF OST400Set01/6

- 1 Intensiv Ortho-Strips Tray REF OST400
- 1 Intensiv Ortho-Strips System, Opener, REF OS08OP-DS
- 2 Intensiv Ortho-Strips System, Central, Medium, REF OS40M-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Coarse, REF OS60C-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Fine, REF OS25F-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Polishing, REF OS15POL-CEN-DS

### EXTENDED

Opening (orange), contoured reduction (red) and fast reduction (brown), finishing (white), polishing (yellow)

### REF OST400Set02/6

- 1 Intensiv Ortho-Strips Tray REF OST400
- 1 Intensiv Ortho-Strips System, Opener, REF OS08OP-DS
- 1 Intensiv Ortho-Strips System, ExtraCoarse, REF OS80XC-DS
- 1 Intensiv Ortho-Strips System, Central, Coarse, REF OS60C-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Medium, REF OS40M-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Fine, REF OS25F-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Polishing, REF OS15POL-CEN-DS

### COMPLETE

Opening (orange), contoured reduction (red) and fast reduction (brown), very fast reduction (blue), finishing (white), polishing (yellow)

### REF OST400Set03/6

- 1 Intensiv Ortho-Strips Tray REF OST400
- 2 Intensiv Ortho-Strips System, Opener, REF OS08OP-DS
- 3 Intensiv Ortho-Strips System, Central, Medium, REF OS40M-CEN-DS
- 1 Intensiv Ortho-Strips System, Central, Polishing, REF OS15POL-CEN-DS

### STANDARD

Opening (orange), contoured reduction (red), polishing (yellow)

## Indications

- Set03: Ortho-Strips standard protocol; opening, reducing/ measuring and polishing
- Set02: Ortho-Strips Complete selection
- Set01: Ortho-Strips standard protocol incl. 60µm (brown) for faster abrasion and 25µm (white) for subsequent finishing

## Benefits

- Arrangement of Intensiv Ortho-Strips during Orthodontic stripping procedures
- Easy introduction into the interdental contact points
- Prevention of cervical step formation
- Dentine hypersensitivity would not be supported