

Intensiv ProxoPolish

Pat. 699 819

Unique one-sided diamond-coated metal strip for efficient stain removal and manual polishing of proximal surfaces

Indications

- Stain removal in proximal areas
- Gloss polishing in proximal areas
- To complete the professional dental hygiene in the proximal zone

Benefits

- Efficient manual cleaning and polishing in one treatment process with the same instrument
- Ergonomic strip handles and short strip for simple and easy use

Dr. A. Devigus, Bülach, Switzerland



Product Description

- Metal strip, diamond-coated on one side only with two working parts of different grit size
- Perforated anatomical ends for a better gripping
- Sterilizable and reusable
- Total length: 80 mm | Height: 2,5 mm



REF PXPO1508/6

■ 15 µm ■ 8 µm