

# Intensiv ProxoStrip Anterior<sup>Pat. 699 819</sup>

The only short metal diamond-coated strip with grips  
for interproximal finishing of anterior restorations

## Indications

- Finishing and polishing of fillings and crown borders in anterior proximal areas
- Removal of excess material and overhangs in filling surfaces and crown borders
- Bevelling in adhesive restorations in the anterior proximal cervical areas

## Benefits

- Finishing and polishing in one process
- Ergonomic grip for simple, safe use
- Short strip for improved finger support while stripping and shorter leverage with reduced distance to the operating area

Dr. A. Devigus, Bülach, Switzerland



## Product Description

- Metal strip, diamond-coated on one side only with two working parts of different grit size
- Perforated anatomical ends for a better gripping
- Sterilizable and reusable
- Total length: 55 mm | Height: 2,5 mm



REF PXA4015/6, REF PXA4015/12    ■ 40 µm    ■ 15 µm