

# Intensiv ProxoStrip Plus Pat. 699 819

A unique, diamond-coated metal strip with grips in plastic for the final polishing of proximal surfaces


## Indications

- Finishing and polishing of fillings and crown borders in proximal areas
- Removal of excess material and overhangs in filling surfaces and crown borders

## Benefits

- Finishing and polishing in one process
- Short strip for improved finger support
- Clinically perfectly polished proximal surfaces

Dr. A. Devigus, Bülach, Switzerland



## Product Description

- Metal strip, diamond-coated on one side only with two working parts of different grit size
- Anatomical ends in sterilizable plastic for a better gripping
- Sterilizable and reusable
- Total length: 83 mm | Height: 2,5 mm



REF PXP4015/6, REF PXP4015/12

■ 40 µm ■ 15 µm