

Intensiv UniglossPaste

Universal single step diamond high gloss polishing paste with micro-grain specifically adapted to the surface

Indications

- Final polishing and high gloss shining for composites

Benefits

- Can be used with either standard nylon bristle junior cup brush or Prophylaxis rubber cup
- High gloss polishing in only 15 seconds
- Same high gloss on restored and natural surfaces
- Only one paste for high gloss polishing

Dr. Simona Giani, Varese, Italy



Product Description

- Polishing paste with a significant proportion of ultra-fine natural diamonds, adapted in its composition to the surface roughness of composites
- The selected mixture of diamond granules corresponds to the the smallest dimensions of the ultrafine grits



REF 2500 - Syringe of 2.5g