

Intensiv UniglossPolisher^{Pat. 2564803}

Shape-stable polymer diamond instrument for high polishing of aesthetic restorations

Indications

- Final high gloss polishing of all accentuated convex and concave surfaces of aesthetic restorations
- Maintenance of existing aesthetic restorations

Benefits


- Stable form of flame and cup shape through multiple use
- Superior gloss effect and minimal surface roughness
- Complete fissure polishing thanks to flexible but stable tip
- Protection of marginal closure of restorations

Dr. Simona Giani, Varese, Italy



Product Description

- Diamond Polymer Polisher highly loaded with ultra-fine diamonds
- Slim flame and compact cup, uncoloured
- Inner core of flame and cup shape with mandrel, consisting of resistant polymer
- Sterilizable, reusable



REF 2300/6 - Package of 3 pcs each
REF 2310/6 - 6 flames | REF 2320/6 - 6 cups